

Data Sheet

SILICATEIN ALPHA (*SUBERITES DOMUNCULA*) ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | AC1011.1 / AC1011.2 |
| Immunogen: | Recombinant silicatein α from <i>Suberites domuncula</i> |
| Swiss-Prot No: | Q2MEV3 |
| Gene Information: | Gene Name: silicaa-g |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Silicatein- α and silicatein- β , tested against following species: <i>Suberites domuncula</i> , <i>Lubomirskia baikalensis</i> , <i>Tethya aurantium</i> , <i>Geodia cydonium</i> |
| Contents: | 20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest. |
| Known applications: | ELISA (1:7000), Western Blot (1:1000), immunohistochemistry (1:100) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided |
| References: | <ol style="list-style-type: none">1. Müller WEG, Krasko A, Le Penec G, Schröder HC (2003). Biochemistry and cell biology of silica formation in sponges. <i>Microsc. Res. Tech.</i> 62: 368-377.2. Müller WEG, Rothenberger M, Boreiko A, Tremel W, Reiber A, Schröder HC (2005). Formation of siliceous spicules in the marine demosponge <i>Suberites domuncula</i>. <i>Cell Tissue Res.</i> 321: 285-297.3. Belikov SI, Kaluzhnaya OV, Schröder HC, Krasko A, Müller IM, Müller WEG (2005). Expression of silicatein in spicules from the Baikalian sponge <i>Lubomirskia baicalensis</i>. <i>Cell Biol Int.</i> 29:943-951. |
| Last updated on: | 18 August 2016 |

For research use only

Publishing research using AC1011? Please let us know so that we can cite your publication as a reference.

